



SURGICAL PRODUCT
CATALOGUE 2018

AUSTRALIA &
NEW ZEALAND



CONTENTS

SURGICAL GOWNS

AERO CHROME*	page 4
MICROCOOL*	page 5
AERO BLUE*	page 6
Surgical Gowns	page 7

SURGICAL DRAPES & PACKS

Cardiac Catheterisation / Radiology Drapes	page 8
Cardiovascular Drapes	page 10
ENT Drapes & Packs	page 12
Equipment Covers	page 15
General Drapes & Packs:	
Accessories	page 16
Basic Packs	page 17
Basic Drapes	page 18
Breast / Chest Drapes	page 20
Fenestrated Drapes	page 21
Laparoscopy / Laparotomy Drapes	page 22
Minor Procedure Drapes	page 24
Thyroid Drapes	page 25
Universal Packs	page 26
Universal Drapes	page 29
Neuro / Spine Drape	page 31
Obstetrics and Gynaecology Drapes & Packs:	
C-Section Drapes	page 32
Laparoscopy Packs	page 33
Laparoscopy Drapes	page 34
Lithotomy Packs	page 35
Lithotomy Drapes	page 36
Under Buttocks Drapes	page 37
Ophthalmic Drapes & Packs	page 38
Orthopaedic Drapes & Packs:	
Ortho Packs	page 39
Extremity Drapes	page 41
Hand Drapes	page 43
Hip Drapes	page 44
Knee Arthroscopy Drapes	page 45
Shoulder Drapes	page 47
U & Bar Drapes	page 49
Urology Drapes & Packs	page 53
Microbial Sealant	page 54
Index	page 55




SURGICAL GOWNS



AERO CHROME* BREATHABLE PERFORMANCE SURGICAL GOWN

- Soft Knitted V-neck Collar
- Clearly Identifiable Silver Colour
- Soft, Light and Breathable Fabric
- Raglan Sleeves
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Blue Huck Towel Included for all Sterile Codes
- AAMI Level 4 / EN 13795 High Performance

Code	Description	Packaging
44680	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Small	Items/Case = 30
44661	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
44681	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large	Items/Case = 30
44662	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
44682	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
44664	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
44683	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
44666	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
44684	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24
44685	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
44667	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
44686	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
44668	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
44687	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
44669	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30

-  Non-Reinforced Fabric
-  Fabric-Reinforced
-  High Performance

All codes in **bold** indicate a sterile product.
NOTE: Please see package inserts for directions for use for each of these products.

SURGICAL GOWNS




MICROCOOL* BREATHABLE HIGH PERFORMANCE SURGICAL GOWNS WITH SECURE-FIT*

- Innovative Patented Glove-gown Interface
- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Blue Huck Towel Included for all Sterile Codes
- AAMI Level 4 / EN 13795 High Performance



Code	Description	Packaging
13359	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, Small	Items/Case = 30
92320	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, Large	Items/Case = 30
72454	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, Large	Items/Case = 40
92321	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, X-Large	Items/Case = 28
30734	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, X-Large	Items/Case = 384
92323	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, XX-Large	Items/Case = 26
72448	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, XX-Large	Items/Case = 32





MICROCOOL* BREATHABLE HIGH PERFORMANCE SPECIALTY SURGICAL GOWNS WITH SECURE-FIT*

- Innovative Patented Glove-gown Interface
- Patented Breathable Fabric
- Low Linting Fabric
- Patented A-Line Design for Seated Procedures
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Blue Huck Towel Included for all Sterile Codes
- AAMI Level 4 / EN 13795 High Performance



Code	Description	Packaging
92330	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, Large	Items/Case = 28
92331	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, X-Large	Items/Case = 26
92333	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, XX-Large	Items/Case = 26

-  Non-Reinforced Fabric
-  Fabric-Reinforced
-  High Performance

All codes in **bold** indicate a sterile product.
NOTE: Please see package inserts for directions for use for each of these products.



SURGICAL GOWNS



AERO BLUE* PERFORMANCE SURGICAL GOWNS


- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Green
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Blue Huck Towel Included for all Sterile Codes
- AAMI Level 3 / EN 13795 High Performance

Code	Description	Packaging
41870	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
43941	AERO BLUE* Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
41871	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
41739	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
42620	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41725	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
41872	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
41740	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
42621	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41727	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
41873	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
41741	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
43940	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
43623	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30
41874	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24

-  Non-Reinforced Fabric
-  Fabric-Reinforced
-  High Performance

All codes in **bold** indicate a sterile product.
NOTE: Please see package inserts for directions for use for each of these products.




SURGICAL GOWNS



SURGICAL GOWNS


- EVOLUTION* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Blue Huck Towel Included for all Sterile Codes
- EN 13795 Standard Performance

Code	Description	Packaging
90003	Surgical Gown, Sterile, Small To Medium	Items/Case = 36
90018	Surgical Gown, Sterile, Large	Items/Case = 36
72978	Surgical Gown, Non-Sterile, Large	Items/Case = 640
90019	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
90048	Surgical Gown, Sterile, X-Large	Items/Case = 32
72988	Surgical Gown, Non-Sterile, X-Large	Items/Case = 640
90049	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
70392	Surgical Gown, Non-Sterile, X- Large, X-Long	Items/Case = 60
90050	Surgical Gown, Sterile, XX-Large	Items/Case = 28
90138	Surgical Gowns Double Pack	Items/Case = 18
	Contents:	
72978	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	

-  Non-Reinforced Fabric
-  Fabric-Reinforced
-  High Performance

All codes in **bold** indicate a sterile product.
NOTE: Please see package inserts for directions for use for each of these products.

CARDIAC CATHETERISATION / RADIOLOGY DRAPES




PACEMAKER DRAPE, 221cm x 343cm

- Two Clear Side Windows 55cm x 342cm
- CONTROL* Plus Fabric Reinforcement
- Conformable Breather Bar 35cm
- Two Adhesive Fenestrations 15cm x 20cm
- Bilateral Pockets 15cm x 81cm

Code	Description	Packaging
89368	Pacemaker Drape, Sterile	Items/Case = 11
79368	Pacemaker Drape, Non-Sterile	Items/Case = 12

BRACHIAL ANGIOGRAPHY DRAPE, 111cm x 78cm


- CONTROL* Plus Fabric Reinforcement
- 7cm Adhesive Circular Fenestration



Code	Description	Packaging
89471	Brachial Angiography Drape, Sterile	Items/Case = 50

RADIAL FEMORAL ANGIOGRAPHY DRAPE, 240cm x 380cm

- Two 63cm x 350cm Windows
- Two 12cm x 7cm Oval Radial Fenestrations
- Two 12cm Round Femoral Fenestrations



Code	Description	Packaging
89583	HALYARD BASICS* Radial Femoral Angiography Drape, Sterile	Items/Case = 12

RADIAL FEMORAL ANGIOGRAPHY DRAPE, 244cm x 381cm

- Two 7cm Adhesive Circular Fenestrations
- Two 4cm x 6cm Oval Fenestrations
- Two 61cm x 381cm Clear Windows
- CONTROL* Fabric Reinforcement

Code	Description	Packaging
44976	Radial/Femoral Angiography Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


CARDIAC CATHETERISATION / RADIOLOGY DRAPES



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOW, 208cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- One 68cm x 315cm Windows
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89709	Extra Long Femoral Angiography Drape with Window, Sterile	Items/Case = 14
79709	Extra Long Femoral Angiography Drape with Window, Non-Sterile	Items/Case = 16



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOWS, 221cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- Two 68cm x 343cm Windows
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89711	Extra Long Femoral Angiography Drape with Dual Windows, Sterile	Items/Case = 16
13296	Extra Long Femoral Angiography Drape with Dual Windows, Non-Sterile	Items/Case = 180



EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220cm x 342cm

- Extra-Long 220cm x 342cm Film Based Drape
- 81cm x 342cm Absorbent Fabric Reinforcement
- 78cm x 342cm Patient Comfort Layer
- Two 8cm Adhesive Circular Fenestrations
- Two 68cm x 342cm Clear Plastic Windows

Code	Description	Packaging
89749	Extra-Long Femoral Angiography Drape, Sterile	Items/Case = 16
79749	Extra-Long Femoral Angiography Drape, Non-Sterile	Items/Case = 144

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

CARDIOVASCULAR DRAPES




187cm x 241cm

PERINEAL DRAPE, 187cm x 241cm

- Split 44cm x 95cm
- CONTROL* Plus Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
76356	Perineal Drape, Non-Sterile	Items/Case = 17




CVARTS* HEAVY-DUTY UNIVERSAL DRAPE SYSTEM

- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89454	CVARTS* Heavy-Duty Universal Drape System, Sterile	Items/Case = 5

Contents:

- 1-CVARTS* Top Drape, 254cm x 152cm
- 1-CVARTS* Bottom Drape, 193cm x 152cm
- 2-CVARTS* Side Drapes, 96cm x 213cm, Longitudinal Flaps with Hook and Loop



CVARTS* HEAVY-DUTY CV SPLIT DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 15cm x 213cm Split
- Longitudinal Flaps with Hook and Loop
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89452	CVARTS* Heavy-Duty CV Split Drape, Sterile	Items/Case = 6
79452	CVARTS* Heavy-Duty CV Split Drape, Non-Sterile	Items/Case = 7

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

CARDIOVASCULAR DRAPES



CVARTS* BILATERAL UNDERLEG / PERINEAL DRAPE, 213cm x 203cm

- Perineal Cover
- Umbilicus Cover
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89453	CVARTS* Bilateral Underleg / Perineal Drape, Sterile	Items/Case = 14
79453	CVARTS* Bilateral Underleg / Perineal Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

ENT DRAPES & PACKS




ENT PACK III

Code	Description	Packaging
88223	ENT Pack III	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89603	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
89431	79430	1-Mayo Stand Cover 80cm x 145cm
89352	79351	1-Suture Bag
		1-Head Drape 91cm x 127cm: <ul style="list-style-type: none"> • Turban • Positioning Tape
		1-ENT Split Drape 193cm x 315cm: <ul style="list-style-type: none"> • with 7 x 71cm Split • Positioning Tape • CONTROL* Fabric Reinforcement • Tube Holders




EENT PACK I

Code	Description	Packaging
88711	EENT Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Absorbent Towel, 38cm x 56cm
		1-Suture Bag
89431	79430	1-Head Drape, 127cm x 91cm
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split



EENT PACK II

Code	Description	Packaging
88716	EENT Pack II	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown with Towel, Large
		1-Microsurgery Mayo Stand Cover
		1-Suture Bag
89431	79430	1-Head Drape, 127cm x 91cm
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split



All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

ENT DRAPES & PACKS



EENT PACK IV


Code	Description	Packaging
88719	EENT Pack IV	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Surgical Gown with Reinforced Sleeves, In Overwrap, with Towel, Large 1-Surgical Gown, with Reinforced Sleeves, Large
89701	79720	1-Absorbent Towel, 38cm x 56cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89438	79438	1-Head Drape, 104cm x 81cm 1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split



EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad


Code	Description	Packaging
89352	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25



HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm

Code	Description	Packaging
89431	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
89438	Head Drape, Sterile	Items/Case = 38

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ENT DRAPES & PACKS



DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm from the top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
89543	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15



EENT FENESTRATED DRAPE, 193cm x 229cm

- 6cm x 5cm Fenestration
- CONTROL Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89421	EENT Fenestrated Drape, Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

EQUIPMENT COVERS



MAYO STAND COVERS

Code	Description	Packaging
88665	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Sterile	Items/Case = 80
78656	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Non-Sterile	Items/Case = 140
89601	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
89602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
14007	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
76358	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 68
89603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44
79603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Non-Sterile	Items/Case = 120



BACK TABLE COVERS

Code	Description	Packaging
88666	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Sterile	Items/Case = 45
78664	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Non-Sterile	Items/Case = 60
89632	Back Table Cover, HALYARD BASICS*, 100cm x 150cm, Sterile	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
42217	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42224	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42306	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
42309	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
42367	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Fan Fold	Items/Case = 160
29258	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20
29259	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 22
29269	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 300
37498	Back Table Cover, Two Tier Heavy Duty, 152cm, Sterile	Items/Case = 10

OVERHEAD / CV TABLE COVERS

Code	Description	Packaging
89621	Overhead Table Cover, 216cm x 193cm, 112cm x 152cm Skirt, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


ACCESSORIES

DISPOSABLE HUCK TOWELS, 51cm x 64cm

Code	Description	Packaging
89021	Disposable Huck Towel, Sterile, Pack of 2	Items/Case = 142
79725	Disposable Huck Towel, Non-Sterile	Items/Case = 300

TAPE STRIPS, 10cm x 55cm


Code	Description	Packaging
89053	Tape Strips, Sterile, Pack of 2	Items/Case = 300



DIATHERMY POUCH, 38cm x 41cm

- Adhesive Backing Strip

Code	Description	Packaging
89945	Diathermy Pouch, Sterile	Items/Case = 150
79945	Diathermy Pouch, Non-Sterile	Items/Case = 220



UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
89731	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS

BASIC PACKS

BASIC PACK VI

Code	Description	Packaging	
88161	Basic Pack VI	Items/Case = 7	
	Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498		1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311		1-ULTRA Surgical Gown, In Overwrap with Towel, Large
95121	71321		1-ULTRA Surgical Gown, with Towel, X-Large
89601	79601		1-Mayo Stand Cover, Reinforced, 58cm x 137cm
89731	79731		4-Utility Drapes, with Tape
			1-Suture Bag
89101	79101		1-Half Drape, 152cm x 112cm
89121	79121		1-Large Drape, 152cm x 193cm

BASIC PACK XV

Code	Description	Packaging	
88178	Basic Pack XV	Items/Case = 15	
	Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498		1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89601	79601		1-Mayo Stand Cover, Reinforced, 58cm x 137cm
			1-Suture Bag
89121	79121		1-Large Drape, 152cm x 193cm

SET-UP PACK I

Code	Description	Packaging	
88171	Set-Up Pack I	Items/Case = 14	
	Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498		1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311		1-Ultra Surgical Gown In Overwrap with Towel, Large
95121	71321		1-ULTRA Surgical Gown, X-Large
89601	79601		1-Mayo Stand Cover, Reinforced, 58cm x 137cm
			1-Suture Bag

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS

BASIC DRAPES




UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
89353	Utility Drape, Sterile	Items/Case = 120
76323	Utility Drape, Non-Sterile	Items/Case = 200




HALF DRAPE, 152cm x 112cm

Code	Description	Packaging
89101	Half Drape, Sterile	Items/Case = 54
79106	Half Drape, Non-Sterile	Items/Case = 1000




MEDIUM DRAPE, 193cm x 112cm

Code	Description	Packaging
89111	Medium Drape, Sterile	Items/Case = 50
79112	Medium Drape, Non-Sterile	Items/Case = 800



MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
89023	Medium Drape with Tape, Sterile	Items/Case = 50



LARGE DRAPE, 152cm x 193cm

Code	Description	Packaging
89121	Large Drape, Sterile	Items/Case = 36
79120	Large Drape, Non-Sterile	Items/Case = 600

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



THREE-QUARTER DRAPE, 152cm x 193cm


- CONTROL* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
89141	Three-Quarter Drape, Sterile	Items/Case = 20



X-LARGE DRAPE, 193cm x 254cm


Code	Description	Packaging
89131	X-Large Drape, Sterile	Items/Case = 20
79132	X-Large Drape, Non-Sterile	Items/Case = 24



SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89828	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89343	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50
79343	Standard Absorbent Drape with Tape, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


BREAST / CHEST DRAPES



CHEST U-DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders

Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22



BREAST DRAPE, 254cm x 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89236	Breast Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


FENESTRATED DRAPES



EPIDURAL DRAPE WITH TAPE, 50cm x 75cm

- 8cm x 10cm Oval Fenestration, Centred
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
89636	HALYARD BASICS* Epidural Drape with Tape, Sterile	Items/Case = 40



VARIABLE FENESTRATED DRAPE, 75cm x 100cm

- Variable Fenestration
- Half Circle On Each Drape with 12cm Diameter
- Spunbond Film
- Positioning Tape

Code	Description	Packaging
89680	HALYARD BASICS* Variable Fenestrated Drape, Sterile	Items/Case = 30



CV LINE DRAPE, 160cm x 240cm

- CONTROL* Fabric Reinforcement
- 12cm x 12cm Fenestration

Code	Description	Packaging
89835	CV Line Drape, Sterile	Items/Case = 10
79835	CV Line Drape, Non-Sterile	Items/Case = 14



CV FULL BODY DRAPE, 173cm x 279cm

- CONTROL* Plus Fabric Reinforcement
- 10cm Circular Adhesive

Code	Description	Packaging
31445	CV Full Body Drape, Non-Sterile	Items/Case = 27

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


LAPAROSCOPY / LAPAROTOMY DRAPES



LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
89229	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers


Code	Description	Packaging
89233	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15
79233	Laparoscopic Abdominal Drape, Non-Sterile	Items/Case = 18



LAPAROSCOPIC DRAPE, 264cm x 89cm x 315cm

- 31cm x 33cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Clear Film Pouch
- Tube Holders
- Armboard Covers
- Under Leg Wings

Code	Description	Packaging
29240	Laparoscopic Drape, Sterile	Items/Case = 15
29241	Laparoscopic Drape, Non-Sterile	Items/Case = 19



HALYARD BASICS* LAPAROSCOPIC/CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 31cm x 33cm Abdominal Fenestration
- Hook and Loop Tube and Line Control
- Fabric Reinforcement
- Armboard Covers

Code	Description	Packaging
89538	Laparoscopic/Cholecystectomy Drape, Sterile	Items/Case = 25
79538	Laparoscopic/Cholecystectomy Drape, Non-Sterile	Items/Case = 26

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89241	Paediatric Laparotomy Drape, Sterile	Items/Case = 15
79240	Paediatric Laparotomy Drape, Non-Sterile	Items/Case = 25



LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89282	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15



ROBOTIC ABDOMINAL / PELVIC DRAPE, 264cm x 406cm

- 30cm x 10cm x 48cm Trapezoidal Abdominal Fenestration
- Clear Film Poly Pouches
- Tube Holders and Line Control
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
29283	Robotic Abdominal/Pelvic Drape, Sterile	Items/Case = 8
29284	Robotic Abdominal/Pelvic Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


MINOR PROCEDURE DRAPES



MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89206	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89207	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89208	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89209	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23



CIRCUMCISION DRAPE, 31cm x 51cm

- 2.5cm Circular Fenestration
- CONTROL* Fabric Reinforcement

Code	Description	Packaging
44972	Circumcision Drape, Sterile	Items/Case = 40

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS

THYROID DRAPES



THYROID DRAPE, 193cm x 315cm

- 12cm x 12cm Diamond Fenestration
- CONTROL* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders

Code	Description	Packaging
89255	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


UNIVERSAL PACKS

UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES



Code	Description	Packaging
88759	Universal Pack – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 15cm x 56cm 1-Suture Bag
34902	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
34906	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
34904	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm

UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES



Code	Description	Packaging
88769	Universal Pack – Minor	Items/Case = 14
Sterile Component Code	Non-Sterile Component Code	Contents
89611	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
	79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
	79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
	79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS



UNIVERSAL PACK I

Code	Description	Packaging
88878	Universal Pack – I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone Reinforced, 112cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm
	79756	2-HALYARD BASICS* Side Drape, 193 x 112cm
	79754	1-HALYARD BASICS* Top Drape, 274 x 127cm
	79755	1-HALYARD BASICS* Bottom Drape, 193 x 173cm



UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging
88774	Universal Pack – III	Items/Case = 8
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS



U-BAR* PACK

Code	Description	Packaging
88334	U-BAR* Pack	Items/Case = 16
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm



U-BAR* PACK – IV

Code	Description	Packaging
88773	U-BAR* Pack – IV	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm

HALYARD BASICS* U-BAR PACK I

Code	Description	Packaging
88875	HALYARD BASICS* U-BAR Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89532	79532	1-HALYARD BASICS* U-Drape, 203cm x 315cm, with 10cm x 112cm Split
89533	79533	1-HALYARD BASICS* Bar Drape, 244cm x 112cm
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 112cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


UNIVERSAL DRAPES



TOP DRAPE, 273cm x 127cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
34906	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



BOTTOM DRAPE, 193cm x 195cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34902	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



SIDE DRAPE, 195cm x 105cm

- CONTROL* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34904	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



TOP DRAPE, 276cm x 150cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769


Code	Description	Packaging
79095	Top Drape, Non-Sterile	Items/Case = 60



BOTTOM DRAPE, 165cm x 193cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79096	Bottom Drape, Non-Sterile	Items/Case = 70



SIDE DRAPE, 114cm x 84cm

- CONTROL* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79097	Side Drape, Non-Sterile	Items/Case = 200

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

NEURO / SPINE DRAPE



CRANIOTOMY DRAPE WITH POUCH, 310cm x 188cm x 340cm

- CONTROL* Plus Fabric Reinforcement
- 20cm x 30cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port

Code	Description	Packaging
89062	Craniotomy Drape with Pouch, Sterile	Items/Case = 9
79075	Craniotomy Drape with Pouch, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS




C-SECTION DRAPES

C-SECTION DRAPE WITH CLEAR SCREEN, 264cm x 193cm x 305cm

- Triangular Fenestration
- 360° Fluid Collection Pouch
- Clear Screen
- Tube Holders
- Armboard Covers


Code	Description	Packaging
44968	C-Section Drape with Clear Screen, Sterile	Items/Case = 6
44967	C-Section Drape with Clear Screen, Non-Sterile	Items/Case = 100



C-SECTION DRAPE WITH POUCH AND FENESTRATION, 264cm x 193cm x 305cm

- 36cm x 31cm Adhesive Area with Fenestration
- 360° Fluid Collection Pouch
- CONTROL* Fabric Reinforcement
- Tube Holders
- Armboard Covers


Code	Description	Packaging
44978	C-Section Drape with Pouch and Fenestration, Sterile	Items/Case = 6
44977	C-Section Drape with Pouch and Fenestration, Non-Sterile	Items/Case = 100



C-SECTION DRAPE WITH FENESTRATED INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area with Triangular Fenestration
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89528	C-Section Drape with Fenestrated Incise, Sterile	Items/Case = 8
79528	C-Section Drape with Fenestrated Incise, Non-Sterile	Items/Case = 150



C-SECTION DRAPE WITH FULL INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89529	C-Section Drape with Full Incise, Sterile	Items/Case = 8
79529	C-Section Drape with Full Incise, Non-Sterile	Items/Case = 150

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

LAPAROSCOPY PACKS

LAPAROSCOPY PACK I



Code	Description	Packaging
88088	Laparoscopy Pack I	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89411		1-Under Buttocks Drape, Plastic, 75cm x 112cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89119	79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, 249cm x 104cm x 251cm <ul style="list-style-type: none"> • 15cm x 20cm Rectangular Abdominal Fenestration • 7cm x 15cm Perineal Fenestration with Flap • Tube Holders



LAPAROSCOPY / PELVISCOPY PACK

Code	Description	Packaging
88222	Laparoscopy/Pelviscopy Pack	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89413		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
89408	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm 1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm <ul style="list-style-type: none"> • 18cm x 23cm Trapezoidal Abdominal Fenestration • 8cm x 15cm Perineal Fenestration with Flap

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LAPAROSCOPY DRAPES



LAPAROSCOPY / PELVISCOPY DRAPE WITH FLUID COLLECTION POUCH, 284cm x 76cm x 314cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Line Controls
- Armboard Covers

Code	Description	Packaging
89219	Laparoscopy/Pelviscopy Drape Drape with Pouch, Sterile	Items/Case = 10



LAPAROSCOPY / PELVISCOPY DRAPE, 254cm x 76cm x 312cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Tube and Line Controls

Code	Description	Packaging
89224	Laparoscopy/Pelviscopy Drape, Sterile	Items/Case = 10
79224	Laparoscopy/Pelviscopy Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

LITHOTOMY PACKS



LITHOTOMY PACK II

Code	Description	Packaging
88521	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89413	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89101	79101	1-Half Drape, 152cm x 111cm



LITHOTOMY PACK V

Code	Description	Packaging
88551	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89701	79720	1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LITHOTOMY DRAPES



LITHOTOMY DRAPE, 254cm x 76cm x 221cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Attached Leggings


Code	Description	Packaging
89203	Lithotomy Drape, Sterile	Items/Case = 16



GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89204	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10



LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55



LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
89408	Leggings, Sterile Pair	Items/Case = 26
79408	Leggings, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


UNDER BUTTOCKS DRAPES



UNDER BUTTOCKS DRAPE, ZONE-REINFORCED, MEDIUM, 96cm x 112cm


- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
79411	Under Buttocks Drape, Zone-Reinforced, Medium, Non-Sterile	Items/Case = 115



UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL* Plus Absorbent Pad
- 72cm x 66cm Triangular Pouch

Code	Description	Packaging
89414	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30
79414	Under Buttocks Drape, with Fluid Collection Pouch, Non-Sterile	Items/Case = 80



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure
- 69cm x 61cm x 7cm Triangular Pouch

Code	Description	Packaging
89415	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


OPHTHALMIC DRAPES & PACKS



BASIC OPHTHALMIC DRAPE, 132cm x 142cm

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
89304	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100



OCUARTS* FULL BODY OPHTHALMIC DRAPE, 183cm x 292cm

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders

Code	Description	Packaging
89356	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS

ORTHOPAEDIC PACKS

ORTHO PACK I – MAJOR



Code	Description	Packaging
88411	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
92321	30734	2-MICROCOOL* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89331	79330	1-Split Drape, 152cm x 178cm with 10cm x 53cm Split 1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split
89381	79360	1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm



ORTHOARTS* HIP PACK WITH POCKETS

Code	Description	Packaging
88452	ORTHOARTS* Hip Pack with Pockets	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, in Overwrap, with Towel, Large
92321	30734	2-MICROCOOL* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89821	79515	1-Stockinette, Large
89121	79121	1-Large Drape, 152cm x 193cm
89333	79333	1-Split Drape, Large, 152cm x 205cm with 10cm x 79cm Split
89345	79345	1-ORTHOARTS* Hip Drape with CONTROL* Plus Reinforcement, 289cm x 223cm x 320cm, with 21cm x 17cm Conformable Teardrop Fenestration 1-Yankauer Suction Tip 1-304cm Suction Tubing

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* LOWER EXTREMITY

Code	Description	Packaging
88471	ORTHOARTS* Lower Extremity Pack	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large 1-Extremity Prep Pad
95221	74921	1-ULTRA Fabric-Reinforced Surgical Gown, with Towel, X-Large
92331		1-MICROCOOL* Breathable High Performance Specialty Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-304cm Suction Tubing
89821	79515	1-Stockinette, Large
89331	79330	1-Split Drape, 152cm x 178cm with 10cm x 53cm Split
89278	13233	1-ORTHOARTS* Lower Extremity Drape with CONTROL* Plus Reinforcement, with 6cm Conformable Circular Fenestration, 3-Drape Tape Strips

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


EXTREMITY DRAPES



LOWER EXTREMITY DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89016	Lower Extremity Drape, Sterile	Items/Case = 12



LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89276	Lower Extremity Drape, Sterile	Items/Case = 14



ORTHOARTS* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89278	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84



ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration


Code	Description	Packaging
89279	Universal Extremity Drape, Sterile	Items/Case = 14
79279	Universal Extremity Drape, Non-Sterile	Items/Case = 150



BILATERAL LIMB DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders


Code	Description	Packaging
89291	Bilateral Limb Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE WITH ARMBOARD COVERS, 295cm x 223cm x 366cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89380	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12



UPPER EXTREMITY DRAPE, 193cm x 305cm

- 3cm Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89271	Upper Extremity Drape, Sterile	Items/Case = 14
79273	Upper Extremity Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS

HAND DRAPES




HAND DRAPE, 289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

77450	Hand Drape, Sterile	Items/Case = 11
--------------	---------------------	-----------------

13227	Hand Drape, Non-Sterile	Items/Case = 108
-------	-------------------------	------------------



ORTHOARTS* HAND DRAPE, 193CM X 381CM

- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

89272	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6
--------------	--------------------------------	----------------

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


HIP DRAPES



ORTHOARTS* HIP DRAPE WITH LEG POCKETS, 289cm x 223cm x 320cm

- Leg Pockets On Both Sides
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89345	ORTHOARTS* Hip Drape with Leg Pockets, Sterile	Items/Case = 6
79345	ORTHOARTS* Hip Drape with Leg Pockets, Non-Sterile	Items/Case = 6



ORTHOARTS* HIP DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89346	ORTHOARTS* Hip Drape, Sterile	Items/Case = 7
79346	ORTHOARTS* Hip Drape, Non-Sterile	Items/Case = 7



HIP DRAPE WITH CLEAR LEG POCKETS, 294cm x 223cm x 345cm

- Clear Leg Pockets
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control

Code	Description	Packaging
89349	Hip Drape with Clear Leg Pockets, Sterile	Items/Case = 9
79349	Hip Drape with Clear Leg Pockets, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS

KNEE ARTHROSCOPY DRAPES


KNEE ARTHROSCOPY DRAPE, 228cm x 315cm



- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 8cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89265	Knee Arthroscopy Drape, Sterile	Items/Case = 10
79265	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 150


KNEE ARTHROSCOPY DRAPE WITH ARMBEARD COVERS, 290cm x 223cm x 315cm



- Armboard Covers
- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89371	Knee Arthroscopy Drape with Armboard Covers, Sterile	Items/Case = 10
79371	Knee Arthroscopy Drape with Armboard Covers, Non-Sterile	Items/Case = 13

ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm



- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7

HALYARD BASICS* KNEE ARTHROSCOPY PACK



Code	Description	Packaging
88886	HALYARD BASICS* Knee Arthroscopy Pack	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88660	78660	1-HALYARD BASICS* Film Split Drape, 152cm x 178cm, with 10cm x 53cm Split
88661	78661	1-HALYARD BASICS* Knee Arthroscopy Drape with Fluid Collection Pouch, 229cm x 315cm
89564	79764	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 152cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 95cm x 112cm Pentagonal Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
36930	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10



HALYARD BASICS* KNEE ARTHROSCOPY DRAPE, 229cm x 315cm

- Fluid Collection Pouch with Suction Port
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
88661	Knee Arthroscopy Drape, Sterile	Items/Case = 20
78661	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


SHOULDER DRAPES



ORTHOARTS* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection Pouch
- 95cm x 97cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
36932	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones
- 90cm x 110cm Trapezoidal Pouch

Code	Description	Packaging
89066	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6
79049	Beach Chair Shoulder Drape without Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 11
79061	Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 44



BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection Pouch
- 103m x 63cm x 2cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
89070	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
79070	Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:



CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK

- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- 55cm x 112cm Trapezoidal Pouch
- Tube Holders


Code	Description	Packaging
89490	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12



ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- 79cm x 149cm x 48cm Trapezoidal Pouch
- Heavy-Duty Fabric


Code	Description	Packaging
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89274	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12



HALYARD BASICS* SHOULDER SPLIT DRAPE, 285cm x 152cm

- 13cm x 34cm Split
- Fabric Reinforcement

Code	Description	Packaging
89559	Shoulder Split Drape, Sterile	Items/Case = 30
79759	Shoulder Split Drape, Non-Sterile	Items/Case = 50



HALYARD BASICS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 158cm

- Fluid Collection Pouch with Suction Port
- Tube Holders

Code	Description	Packaging
89558	Shoulder Arthroscopy Drape, Sterile	Items/Case = 30
79758	Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


U & BAR DRAPES



U-DRAPE, 193cm x 305cm

- 10cm x 101cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89301	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190



LARGE U-DRAPE, 223cm x 317cm

- 15cm x 137cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89311	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* U-DRAPE, 223CM X 274CM

- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89347	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10
13294	ORTHOARTS* U-Drape, Non-Sterile	Items/Case = 132

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* BODY SPLIT DRAPE, 223CM X 294CM

- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89275	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12



SPLIT DRAPE, 152cm x 178cm

- Polyethylene Material
- 10cm x 53cm Split


Code	Description	Packaging
89331	Split Drape, Sterile	Items/Case = 42
79330	Split Drape, Non-Sterile	Items/Case = 115



LARGE SPLIT DRAPE, 152cm x 205cm

- Polyethylene Material
- 11cm x 78cm Split

Code	Description	Packaging
89333	Large Split Drape, Sterile	Items/Case = 34
79333	Large Split Drape, Non-Sterile	Items/Case = 100



HALYARD BASICS* BODY SPLIT DRAPE, 224cm x 295cm

- 10cm x 67cm Split
- Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89557	Body Split Drape, Sterile	Items/Case = 25
79757	Body Split Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



ORTHO SPLIT "T" DRAPE, 297cm x 210cm x 342cm

- Armboard Covers
- 10cm x 53cm Split
- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
36931	Ortho Split "T" Drape, Sterile	Items/Case = 10



ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89381	Ortho Bar Drape, Sterile	Items/Case = 18
79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28



VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- 42cm x 84cm x 17cm Trapezoidal Pouch
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
77723	Vertical Isolation Drape, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS



ISOLATION SCREEN, 170cm x 290cm

- Clear Plastic Drape 170cm x 290cm
- Adhesive Backing Strips

Code	Description	Packaging
89093	Isolation Screen, Sterile	Items/Case = 24
79093	Isolation Screen, Non-Sterile	Items/Case = 30



STOCKINETTES

Code	Description	Packaging
89811	Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
79505	Stockinette, Non-Sterile, Medium, 23cm x 112cm	Items/Case = 75
89821	Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
79515	Stockinette, Non-Sterile, Large, 30cm x 122cm	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


UROLOGY DRAPES & PACKS



238 cm x 175 cm x 170 cm


T.U.R. PACK

Code	Description	Packaging
77885	T.U.R. Pack	Items/Case = 12
	Non-Sterile Component Code	Contents
		1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Absorbent Towel, 38cm x 56cm
76885	1-T.U.R. Drape, 289cm x 175cm x 173cm:	
		<ul style="list-style-type: none"> • Attached Leggings • Fluid Collection Pouch with Screen, Suction Port and Tie • 54cm x 75cm Triangular Pouch • Finger Cot • 7cm Fenestration • 25cm x 20cm Fenestration



CYSTOSCOPY PACK I

Code	Description	Packaging
88611	Cystoscopy Pack I	Items/Case = 12
	Sterile Component Code	Non-Sterile Component Code
		Contents
		1-Back Table Cover, Zone-Reinforced, 112cm x 198cm
		1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuff, 78cm x 122cm
		1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm with 7cm x 15cm Perineal Fenestration and 23cm x 23cm Screen



PERCUTANEOUS NEPHROLITHOTOMY DRAPE (PCNL), 299cm x 182cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Fluid Collection Pouch
- 67cm x 61cm x 7cm Trapezoidal Pouch
- 18cm x 28cm Fenestration

Code	Description	Packaging
42427	Percutaneous Nephrolithotomy Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

MICROBIAL SEALANT



HALYARD* INTEGUSEAL* MICROBIAL SEALANT

- Ergonomic Single-Use Applicator
- Evenly Covers Incision Area in a Single Stroke
- Widely-Used Cyanoacrylate Forms a Safe, Fast-Drying, Breathable Barrier
- Compatible with Most Common Surgical Techniques and Materials
- Need Not be Removed for Skin Closure
- Can be Used with a Variety of Preps and Incise Drapes
- Nonflammable, Latex-Free, Non-Irritating and Non-Cytotoxic

Code	Description	Coverage Area (approx)	Example of Procedures	Packaging
33736	IS 100	10" x 10" 25cm x 25cm	Cranial and Spinal Procedures, Laparotomies, Hernia Repair, Thyroid Surgery, Cardiac Catheterisation	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20
33737	IS 200	20" x 10" 50cm x 25cm	Cholecystectomy, Open Heart, Lung Surgery, Total Joint Replacement	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20

INDEX

Code	Page	Code	Page	Code	Page	Code	Page	Code	Page
13227	43	43941	4	79273	42	88611	53	89347	49
13231	49	44661	4	79275	50	88661	46	89349	44
13233	41	44662	4	79279	42	88665	15	89352	13
13294	49	44664	4	79282	23	88666	15	89353	18
13296	9	44666	4	79304	38	88711	12	89356	38
13359	5	44667	4	79330	50	88716	12	89368	8
14007	15	44668	4	79333	50	88719	13	89371	45
29219	41	44669	4	79335	20	88759	26	89380	42
29220	41	44680	4	79343	19	88769	26	89381	51
29240	22	44681	4	79345	44	88774	27	89408	36
29241	22	44682	4	79346	44	88773	28	89413	37
29283	23	44683	4	79349	44	88875	28	89414	37
29284	23	44684	4	79351	12	88878	27	89415	37
29258	15	44685	4	79356	38	88886	45	89421	14
29259	15	44686	4	79360	51	89016	41	89431	13
29269	15	44687	4	79363	49	89021	16	89438	13
30734	5	44967	32	79368	8	89023	18	89452	10
31445	21	44968	32	79371	45	89053	16	89453	11
33736	54	44972	24	79408	36	89062	31	89454	10
33737	54	44976	8	79411	37	89066	47	89471	8
34901	29	44977	32	79414	37	89070	47	89490	48
34902	29	44978	32	79415	37	89093	52	89492	48
34903	29	70392	7	79430	13	89101	18	89528	32
34904	29	72448	5	79452	10	89111	18	89529	32
34905	29	72454	5	79453	11	89121	18	89538	22
34906	29	72978	7	79482	15	89131	19	89543	14
36930	46	72988	7	79490	48	89141	19	89557	50
36931	51	76323	18	79492	48	89203	36	89558	48
36932	47	76356	10	79505	52	89204	36	89559	48
37498	15	76358	15	79515	52	89206	24	89583	8
41725	5	77450	43	79528	32	89207	24	89601	15
41727	5	77723	51	79529	32	89208	24	89602	15
41739	5	77885	53	79538	22	89209	24	89603	15
41740	5	78656	15	79551	36	89219	34	89621	15
41741	5	78661	46	79602	15	89224	34	89632	15
41870	6	78664	15	79603	15	89229	22	89636	21
41871	5	79049	47	79607	15	89233	22	89680	21
41872	5	79061	47	79709	9	89236	20	89709	9
41873	5	79070	47	79725	16	89241	23	89711	9
41874	5	79075	31	79731	16	89255	25	89731	16
42217	15	79093	52	79749	9	89265	45	89749	9
42218	15	79095	30	79757	50	89267	45	89811	52
42224	15	79096	30	79758	48	89271	42	89821	52
42225	15	79097	30	79759	48	89272	43	89828	19
42231	15	79106	18	79828	19	89274	48	89835	21
42233	15	79112	18	79835	21	89275	50	89945	16
42306	15	79120	18	79945	16	89276	41	90003	7
42307	15	79132	19	88088	33	89278	41	90018	7
42309	15	79204	36	88161	17	89279	42	90019	7
42310	15	79206	24	88171	17	89282	23	90048	7
42366	15	79207	24	88178	17	89291	45	90049	7
42367	15	79208	24	88222	33	89301	49	90050	7
42426	47	79209	24	88223	12	89304	38	90138	7
42427	53	79224	34	88334	28	89311	49	92320	5
42429	45	79233	22	88411	39	89331	50	92321	5
42620	5	79240	23	88452	39	89333	50	92323	5
42621	5	79255	25	88471	40	89343	19	92330	5
43623	5	79265	45	88521	35	89345	44	92331	5
43940	5	79267	45	88551	35	89346	44	92333	5



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



Halyard Australia Pty. Limited
ABN No 25 100 844 488
Level 7, 52 Alfred Street, Milsons Point NSW
2061
Phone: 1800 664 227
Email: CustomerService_ANZ@hyh.com

www.halyardhealth.com.au



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2018 HYH. All rights reserved.
SS:1803:WFHY770