



PATIENT CARE
ACCESSORIES
CATALOGUE

AUSTRALIA &
NEW ZEALAND



QUICK LINKS

COLD THERAPY

Ice Packs

HALYARD* Bilateral Facial Ice Pack	page 96
HALYARD* EyeCare (EENT) Ice Pack	page 96
HALYARD* STAY-DRY Ice Pack	page 96

I.V. SUPPORT

I.V. Arterial And Wrist Supports

HALYARD* Neonatal & Infant/Child I.V. Support	page 97
HALYARD* HAND-AID* Arterial Wrist Support	page 97
HALYARD* HAND-AID* Wrist Support	page 97
HALYARD* Antecubital I.V. Support	page 97

O.R. ACCESSORIES

Support Straps

HALYARD* Stirrup Strap	page 98
HALYARD* TS-30* Personal Utility Strap	page 98
HALYARD* Tourniquet Cover	page 98

PATIENT LINEN

Single Use Linen

HALYARD* Stretcher Sheet	page 99
HALYARD* Draw Sheet	page 99
HALYARD* Pillowcase	page 99

COLD THERAPY

ICE PACKS



HALYARD* BILATERAL FACIAL ICE PACK

- 2 ice packs bonded together to evenly deliver cold therapy to both sides of the jaw
- Can be used for recovery after maxillofacial and dental procedures
- Single-patient use

Code	Description	Size	Packaging
33101	Bilateral Facial Ice Pack	13cm x 30cm	12 eaches/box 2 boxes/case 24 eaches/case



HALYARD* EYECARE (EENT) ICE PACK

- Designed for easy application after plastic surgery, ophthalmic and sinus procedures
- Adjustable strap ensures the pack sits comfortably and securely on the face
- Single-patient use

Code	Description	Size	Packaging
33150	EyeCare (EENT) Ice Pack, 1 Strap	11.4cm x 25cm	15 eaches/box 2 boxes/case 30 eaches/case



HALYARD* STAY-DRY ICE PACK

- For general use
- Delivers up to 2.5 hours of consistent, dry cold therapy solutions
- STAY-DRY construction features 3 layers which help to reduce wetness due to condensation
- 4 ties for versatile application on shoulders, knees, elbows, etc
- Single-patient use

Code	Description	Size	Packaging
33495	STAY-DRY Ice Pack, Small, 4 Ties	13cm x 30cm	25 eaches/box 2 boxes/case 50 eaches/case
33500	STAY-DRY Ice Pack, Large, 4 Ties	16.5cm x 36cm	25 eaches/box 2 boxes/case 50 eaches/case

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

I.V. SUPPORT

I.V. ARTERIAL AND WRIST SUPPORTS



HALYARD* NEONATAL & INFANT/CHILD I.V. SUPPORT

- Padded for comfort
- Aluminium core easily bends to conform to various angles of small wrists and feet
- Single-patient use

Code	Description	Size	Packaging
29965	Neonatal I.V. Support	9cm x 4cm, 13cm straps	30 eaches/case
29967	Infant/Child I.V. Support	16cm x 5cm, 23cm straps	30 eaches/case



HALYARD* HAND-AID* ARTERIAL WRIST SUPPORT

- Uniquely anatomically designed to position the wrist for accurate and efficient application of arterial lines
- The instant support of the foam block saves setup time
- Single-patient use

Code	Description	Size	Packaging
29990	Hand-Aid* Arterial Wrist Support, Paediatric	16cm x 7cm, 23cm straps	20 eaches/case
29980	Hand-Aid* Arterial Wrist Support, Adult	23cm x 9cm, 37cm straps	20 eaches/case
29985	Hand-Aid* Arterial Wrist Support, Replacement Straps & Starter Block, Compatible with 29980		50 eaches/case



HALYARD* HAND-AID* WRIST SUPPORT

- Anatomically designed to position the wrist for accurate and efficient cannulisation in dorsal hand veins
- Single-patient use

Code	Description	Size	Packaging
29970	Hand-Aid* I.V. Wrist Support, Paediatric	14cm x 5cm, 23cm straps	20 eaches/case
29900	Hand-Aid* I.V. Wrist Support, Adult	16cm x 7cm, 23cm straps	20 eaches/case



HALYARD* ANTECUBITAL I.V. SUPPORT

- The unique and comfortable way to immobilise and support during antecubital I.V. therapy
- Single-patient use

Code	Description	Size	Packaging
29920	Antecubital I.V. Support	32cm x 11cm, 48 cm straps	10 eaches/case

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

O.R. ACCESSORIES

SUPPORT STRAPS



HALYARD* STIRRUP STRAP

- Designed to ensure patient comfort in the lithotomy position
- Made of soft fabric and featuring a layer of padding
- Comes with slip ring for easy adjustment
- Single-patient use

Code	Description	Size	Packaging
52712	Personal Stirrup Strap with Slip Ring	48cm x 9cm	10 eaches/pack 4 packs/case 40 eaches/case



HALYARD* TS-30* PERSONAL UTILITY STRAP

- Designed to secure and position the patient during a surgical procedure
- Hook and loop closure for easy adjustment of the strap
- Soft fabric for patient comfort
- Single-patient use

Code	Description	Size	Packaging
55299	TS-30* Personal Utility Strap	76cm x 8cm	12 eaches/pack 2 packs/case 24 eaches/case



HALYARD* TOURNIQUET COVER

- Designed to cover a reusable, pneumatic tourniquet cuff to protect it from skin prep and other procedural fluids
- Soft padding helps prevent pinching and wrinkling of patients' skin
- Comes in two sizes – small and large
- Single-patient use

Code	Description	Size	Packaging
53600	O.R. Tourniquet Cover, Small	53.5cm x 37cm	15 eaches/pack 2 packs/case 30 eaches/case
53605	O.R. Tourniquet Cover, Large	78cm x 37cm	15 eaches/pack 2 packs/case 30 eaches/case

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

PATIENT LINEN



HALYARD* SINGLE USE LINEN

Code	Description	Size	Packaging
67783	Stretcher Sheet, Blue	101.5cm x 213cm (40" x 84"), Bulk	80 eaches/case
67843	Draw Sheet, Poly-Coated*, Blue <small>*Coated with Polyethylene Film</small>	101.5cm x 183cm (40" x 72"), Bulk	80 eaches/case
67803	Pillowcase, Blue	50.5cm x 73.5cm (20in x 29in), Bulk	200 eaches/case

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.